

Premium Dermal Filler
for Facial & Body Rejuvenation

Ultra



Distributor Biowellspharm Co.,Ltd.
Address 3F, 499, Gangseo-ro, Gangseo-gu, Seoul, 07523, Korea
Tel +82-2-3663-2257 / Fax +82-2-3663-2254
E-mail biowp2019@gmail.com / Website www.biowellspharm.co.kr

REGENESS ULTRA FILLER

Product Overview

Regeness Ultra Filler is a premium hyaluronic acid(HA) dermal filler, crafted to restore volume, enhance contours, and reduce facial wrinkles. Each 3mL syringe provides visible results through minimally invasive procedures, leaving the skin smoother and revitalized.

For patients seeking non-surgical breast volume enhancement, Regeness Ultra Filler offers a natural, subtle boost, creating a fuller, more contoured breast line with improved firmness.

Indications for Use

Regeness Ultra Filler is ideal for a range of cosmetic treatments, including :

- **Facial Volume Restoration** : Replenishes lost volume for a more youthful look.
- **Lip Enhancement** : Adds fullness and definition to the lips.
- **Nasolabial Folds Correction** : Smooths out smile lines.
- **Marionette Lines Improvement** : Softens lines extending from the corners of the mouth.
- **Jawline Contouring** : Defines and sculpts the jawline.
- **Non-Surgical Rhinoplasty** : Refines the shape of the nose without surgery.
- **Hand Rejuvenation** : Restores a youthful appearance to the hands.
- **Scar Correction** : Minimizes the appearance of scars for smoother skin.
- **Body Volume Enhancement** : Ideal for contouring and volumizing small areas of the body.

Composition and Key Properties

- **Main Ingredients**

Contains cross-linked hyaluronic acid (24mg/mL) for long-lasting volumizing effects and lidocaine hydrochloride (3mg/mL) to ensure a more comfortable treatment experience.

- **Viscoelasticity**

With a viscoelasticity range of 250-310 Pa, the filler offers a smooth consistency that integrates seamlessly with the skin.

- **Administration**

Regeness Ultra Filler should only be administered by a trained healthcare professional via subcutaneous injection. This ensures optimal outcomes and patient safety.

Post-Treatment Care

Mild side effects such as swelling, redness, or bruising may occur at the injection area.

Patients are advised to avoid extreme temperatures and refrain from applying pressure to treated areas for at least 24 hours post-procedure.

How Supplied

Each set includes three 3mL syringes, providing a total of 9mL of product per package.

Storage and Handling

- Store between 2°C and 25°C, away from direct sunlight.
- Do not freeze, as this may compromise the product's quality.

Sales and Distribution

Regeness Ultra Filler is available for Original Design Manufacturing(ODM) and Original Equipment Manufacturing(OEM) to meet specific client requirements.